


BELLMAN VISIT/RING VIDEO DOORBELL PACK



SKU: WM-Ringkit

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit/Ring Video Doorbell Pack** combines the accessibility of the Bellman Visit alert system with the smart functionality of the Ring Pro Video Doorbell. This integrated package ensures you never miss a visitor, whether you are at home or away. When the Ring Chime sounds, the Bellman Visit Door Transmitter activates your Bellman receivers, alerting you through flashing lights, vibrations, or audible signals — ideal for individuals who are deaf, hard of hearing, or need visual/tactile alerts.

FEATURES

- Combines the Bellman Visit system with the Ring Pro Video Doorbell and Chime for complete home awareness
- Door Transmitter detects the Ring Chime and triggers Bellman Visit receivers for visual, sound, or vibration alerts
- Receive instant mobile notifications and live two-way communication via the Ring app
- 1080p HD wide-angle video with infrared night vision and Colour Pre-Roll for pre-alert footage
- Expandable with other Bellman Visit transmitters and receivers including smoke alarms, alarm clocks, and pagers
- Simple installation with wireless communication between devices

SPECIFICATIONS

- **Included Components:** Bellman Visit Door Transmitter (BE1411), Ring Pro Video Doorbell, Ring Chime
- **Video Quality:** 1080p HD wide-angle with night vision

- **Connectivity:** Wi-Fi (Ring) and radio link (Bellman Visit system)
- **Alerts:** Visual flashes, vibrations, or sound via Bellman Visit receivers
- **Compatibility:** Works with Bellman Visit Flash, Alarm Clock, Pager, and Wrist Receivers

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au