

BELLMAN VISIT SMOKE ALARM PACK



SKU: BE8067

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Bellman Visit Smoke Alarm Pack** provides complete fire safety alerts for individuals who are deaf, hard of hearing, or heavy sleepers. This system offers both visual and tactile notifications through a bright strobe flash and a powerful vibrating bed shaker when smoke is detected. It ensures early warning and reliable communication across the home, allowing users to stay safe and independent at all times.

FEATURES

- Comprehensive smoke alarm alert system with visual and vibrating notifications
- Photoelectric smoke alarm detects smoke and smouldering fires early
- Flash receiver provides bright strobe alerts beside the bed
- Bed shaker delivers strong vibration alerts under the pillow for deep sleepers
- Wireless communication between smoke alarm and receiver for easy setup
- Rechargeable battery backup for continued operation during power outages
- Expandable system compatible with other Bellman Visit transmitters, including doorbell, phone, and baby monitors

SPECIFICATIONS

- **Included Components:** BE1442 Visit Flash (strobe light), BE1289 Visit Smoke Alarm, BE1270 Visit Bed Shaker
- **Alert Type:** Strobe flash and vibration
- **Power Supply:** Flash receiver – mains powered with rechargeable backup; Smoke alarm – sealed 10-year lithium battery

- **Wireless Range:** Typical coverage throughout a standard home environment
- **Compatibility:** Works with additional Bellman Visit transmitters and receivers

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au