



## LIGHT DRIVE PLUS 2.1 POWER ASSIST ADD-ON



**SKU:** 81MM21

## PRODUCT DESCRIPTION

### Purpose & Benefits

The Benoit Systemes Light Drive Plus 2.1 Power Add-On is the heavy-duty version of the popular Light Drive 2.1, engineered in France for users who require extra power, torque, and weight capacity. Fitted with four 130W motors that apply a friction drive to your wheelchair's tyres, it transforms a manual chair into a fully powered system—controlled easily with one hand via a joystick.

With a top speed of 9.98 km/h and a travel range of up to 18 km on a single charge, this add-on is ideal for tackling rougher ground, hills, and longer journeys. Designed to support users up to 240 kg, it attaches and detaches in seconds without tools. Despite its robust performance, the unit remains portable at just 16 kg, and is available in a choice of 15 stylish colours to match your wheelchair.

### Features

- Converts a manual wheelchair into a joystick-controlled powered model
- Four powerful 130W friction drive motors
- Designed to tackle hills, grass, and uneven terrain
- One-handed joystick operation
- Quick, tool-free attachment and removal

- Compatible with most rigid and folding wheelchairs
- Freewheel mode for manual use
- Designed and manufactured in France
- Choice of 15 colour options

## **Specifications**

- Maximum Speed: 9.98 km/h
- Range: Up to 18.02 km
- Unit Weight: 16 kg
- Maximum User Weight: 240 kg
- Motor: 4 × 130W friction drive motors
- Battery: 18 Ah lithium-ion (airline compliant)
- Control Type: Joystick with adjustable driving modes
- Mounting System: Tool-free quick attachment
- Compatibility: Most adult manual wheelchairs
- Power Add-On Type: Rear-mounted friction drive system