


## BETTERLIVING ELEPHANT FEET BED RAISERS & CHAIR RAISERS - SET OF 4

**\$63.00 - \$173.00** Price range: \$63.00 through \$173.00



**SKU:** H-BL16

### PRODUCT DESCRIPTION

#### PURPOSE & BENEFITS

These plastic Elephant Feet Bed Raisers & Chair Raisers are designed to go under the legs of your bed, lounge chair, couch or table to increase the overall height. They allow you to ensure the height of your furniture is right for your needs. Bed Raisers can also be used to raise one end of your bed, to create a slight incline, which can be helpful in managing some digestive, respiratory and circulatory problems.

Made from durable, high impact plastic, these stackable Bed Raisers are available in three different heights. The plastic raiser has a top well that can accommodate square or round legs up to 65 mm wide. They are sold in sets of four. Each set of four can hold up to 250 kilograms of weight, including the furniture and any people.

Please note that your furniture must be sufficiently heavy to remain stable on top of the raisers, so these products are not suitable for lightweight tables or [chairs](#).

The 40mm high raiser has a base width of 125mm. The 100mm high raiser has a base width of 198mm. Combine these two to create a 140mm high raiser with a base width of 198mm.

#### KEY FEATURES

- Durable plastic furniture raisers
- Place under furniture legs to increase the overall height
- Use for beds, couches, tables or chairs
- Suits square and round legs
- Sold in sets of 4

- Holds up to 250kg total

You can buy these elephant legs for beds online using AfterPay and NDIS payment methods.