



## **BETTERLIVING IV POLE - SUITS COMMUNITY BED**

**\$250.00**



**SKU: BE8003**

## **PRODUCT DESCRIPTION**

### **Features and Benefits**

The BetterLiving IV Pole is a sturdy and practical accessory designed to integrate seamlessly with the BetterLiving Community Bed. Ideal for clinical and homecare environments, it ensures essential IV therapy is safely and conveniently accessible. Its height-adjustable design accommodates various patient and treatment needs, while the robust steel construction supports reliable, long-term use.

This accessory attaches directly to the bed frame, maintaining a neat and secure setup without taking up floor space — perfect for improving bedside efficiency and safety in aged care or hospital settings.

### **Key Features**

- Specifically designed for the BetterLiving Community Bed
- Provides safe and secure IV pole attachment to the bed frame
- Height-adjustable for various treatment positions
- Sturdy powder-coated steel construction
- Easy to install and remove for cleaning or repositioning

## Specifications

- **Height Range:** 1100 mm – 1800 mm (adjustable)
- **Mount Type:** Bed-frame attachment