



BETTERLIVING SELF HELP POLE, SUITS COMMUNITY BED AND AVALON ROTATING BED



SKU: BE7003

PRODUCT DESCRIPTION

Purpose & Benefits

This self help pole is designed to increase patient independence by providing a stable support handle for repositioning or sitting up in bed. It's ideal for those with reduced mobility or recovering from surgery, and is compatible with both the BetterLiving Community Bed and the Avalon Rotating Bed.

The pole features a sturdy, powder-coated steel frame with an adjustable nylon strap and ergonomically shaped triangle handle. It attaches securely to the bed frame, giving users a safe, easy-to-grip assist for movement without needing caregiver help.

Features

- Compatible with BetterLiving Community Bed and Avalon Rotating Bed
- Supports safe, independent repositioning
- Triangular ergonomic handle
- Adjustable strap for user comfort
- Durable powder-coated steel construction
- Secure, frame-mounted attachment



• Ideal for aged care, hospitals, and home care

Specifications

• Height: 1400mm

• Product weight: 4.5kg

• Handle type: Triangular ergonomic grip

• Attachment: Bed frame mount

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au