




## BETTERLIVING TOILET SURROUND MAT



**SKU:** RM-7466-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The BetterLiving Toilet Surround Mat provides extended coverage around your toilet for a soft, warm, and secure underfoot surface, reducing the risk of falls in the bathroom.

Its ultra-low profile design fits most toilets and remains in place when using a mobile shower commode. Constructed from fade- and bleach-resistant solution-dyed nylon with a crush-resistant pile, it traps dirt and maintains its performance over time.

The skid-resistant backing repels oil and grease on tiled, timber, or lino surfaces, while an anti-trip border ensures the mat lies flat. Easy to clean by sweeping, vacuuming, or machine washing.

### FEATURES

- Extended coverage for a secure, soft underfoot surface
- Ultra-low profile fits most toilets and mobile shower commodes
- Solution-dyed nylon: fade & bleach-resistant with crush-resistant pile
- Skid-resistant backing repels oil and grease
- Anti-trip border ensures the mat lies flat
- Easy to clean: sweep, vacuum, machine wash, or tumble dry

### SPECIFICATIONS

Length	Width	Opening Width
74 cm	66 cm	30 cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)