

## BISHOP COMFORT BARIATRIC MID BACK LUMBAR SUPPORT CHAIR



**SKU:** AGE500B

### PRODUCT DESCRIPTION

#### PURPOSE & BENEFITS

A fully upholstered, multi-purpose, support chair is designed for use in any environment where an adjustable floor to seat height is required: Lounge, dining, bedroom.

Available in a variety of colour options in both vinyl and fabric.

**Please note:** The standard weight capacity is **350kg**. However, these chairs **can be customised to support heavier Safe Working Loads (SWL)** if required. This becomes a customisation option via **Age Prep Care**.

#### FEATURES

- Bishop Comfort Bariatric Healthcare Chair
- Integrated Lumbar Support
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Floor to Seat Adjustment
- 350kg Weight Capacity

- 10-Year Factory Warranty

## COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

## DIMENSIONS

### Bariatric

Width	650mm
Depth	450mm
Floor to Seat Adjustment	470mm-595mm
Back Rest Height	490mm
Arm Rest Height	170mm
Weight Capacity	350kg
Package Weight	15kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)