



BODYPOINT LEG HARNESS

~~\$264.00~~ Original price was: \$264.00. \$237.60 Current price is: \$237.60.



SKU: WCA8604-H

PRODUCT DESCRIPTION

BENEFITS

Inspired by a therapist's desire for a 5-point harness, the Leg Harness is the go-to solution for people that tend to slide out of their [chairs](#), thrust or have spasticity. It wraps each leg from the rear of the seat to the inner thigh at the front of the chair. This increases the surface area for greater control of movement. The harness has a unique, user-friendly design so that correct use is almost foolproof and it can be used with or without a pelvic belt depending on the specific need. This is a great alternative to a groin strap. Sold as a pair.

KEY FEATURES

- Holds each leg at the inner thigh providing efficient resistance to pelvic rotation and obliquity
- No pressure on the bladder
- Convenient pull straps for adjustment

SPECIFICATIONS

- Small: 230mm
- Medium: 330mm
- Large: 380mm

Buy this BodyPoint Leg Harness online at Medimart Australia, we have NDIS Payment methods as well as AfterPay and ZipPay.