



BREEZY BASIX 2 SELF PROPELLED WHEELCHAIR



SKU: MWSSF

PRODUCT DESCRIPTION

BREEZY BASIX 2 LIGHTWEIGHT WHEELCHAIR

Adjustability and rigidity with a modern look.

The Breezy Basix 2 is a lightweight <u>folding wheelchair</u> that features amazing adjustability and rigidity, combined with a modern and clean look.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest that is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Highly efficient and easy to use - The BasiX 2 is the ideal wheelchair for fleet management.

THIS BREEZY WHEELCHAIR IS;

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion.
- Swing Away and Adjustable Anti-Tips
- Safe working Load 125kg
- Chair weight 14.9kg
- Fixed Backrest

Buy this Breezy Wheelchair online at Medimart using <u>NDIS Payment</u> method and Afterpay, for your shopping convenience.



SPECIFICATIONS

Seat Width: 40.6cm (16"); 45.72cm (18"); 50.8cm (20");

• Seat Depth: 43.5 cm (41 - 46 cm adjustable)

• Seat Height: 42.5 - 51 cm

• Seat Angle: 0º to 6º

• Seat Tilt: no

• Backrest Height: 43.5 cm (41 - 46 cm adjustable)

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au