




BREEZY BASIX 2 SELF PROPELLED WHEELCHAIR

\$922.00 - \$1,123.00 Price range: \$922.00 through \$1,123.00



SKU: MWSSP

PRODUCT DESCRIPTION

BREEZY BASIX 2 LIGHTWEIGHT WHEELCHAIR

Adjustability and rigidity with a modern look.

The Breezy Basix 2 is a lightweight [folding wheelchair](#) that features amazing adjustability and rigidity, combined with a modern and clean look.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest that is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Highly efficient and easy to use - The BasiX 2 is the ideal wheelchair for fleet management.

THIS BREEZY WHEELCHAIR IS;

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage
- Comes with free standard basic foam cushion.
- Swing Away and Adjustable Anti-Tips
- Safe working Load 125kg
- Chair weight 14.9kg
- Fixed Backrest

Buy this Breezy Wheelchair online at Medimart using [NDIS Payment](#) method and Afterpay, for your shopping convenience.

SPECIFICATIONS

- Seat Width: 40.6cm (16"); 45.72cm (18"); 50.8cm (20");
- Seat Depth: 43.5 cm (41 - 46 cm adjustable)
- Seat Height: 42.5 - 51 cm
- Seat Angle: 0° to 6°
- Seat Tilt: no
- Backrest Height: 43.5 cm (41 - 46 cm adjustable)