




BREEZY RUBIX2 XL HEAVY DUTY WHEELCHAIR - BLUE

\$2,170.00



SKU: 074000-001-52-H

PRODUCT DESCRIPTION

Purpose & Benefits of Breezy RubiX² XL Heavy Duty Wheelchair

The **Breezy RubiX² XL Heavy Duty Wheelchair** is a reinforced manual wheelchair designed to provide reliable support, comfort, and adjustability for users requiring a higher weight capacity. Engineered with durability in mind, the RubiX² XL accommodates users up to 170 kg, making it ideal for everyday mobility in both home and care settings. **Please note: Only available in blue.**

With adjustable seat depth, backrest angle, and footplate height, this chair delivers personalised comfort and promotes proper posture. Its heavy-duty frame ensures strength without compromising manoeuvrability, allowing users to move with confidence and independence.

Key Features:

- Heavy-duty frame with a 170 kg weight capacity
- Adjustable seat depth, back angle, and footrest height for tailored comfort
- Reinforced cross-brace and durable construction for long-term reliability
- Swing-away, height-adjustable leg rests and removable armrests
- Large rear wheels and ergonomic hand rims for easy self-propulsion

- Compatible with a range of seating and positioning accessories
- Ideal for home use, rehab centres, and long-term care facilities

Specifications:

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight
074000-001-52	Breezy Rubix 2 XL - 52cm Wide	410 - 510 mm	52cm (20.4")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg
074000-001-56	Breezy Rubix 2 XL - 56cm Wide	410 - 510 mm	56cm (22")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg
074000-001-60	Breezy Rubix 2 XL - 62cm Wide	410 - 510 mm	62cm (24.4")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg