




## **BREEZY RUBIX2 XL HEAVY DUTY WHEELCHAIR**

**\$2,170.00**



**SKU: 074000-001-52-H**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Breezy RubiX<sup>2</sup> XL Heavy Duty Wheelchair**

The **Breezy RubiX<sup>2</sup> XL Heavy Duty Wheelchair** is a reinforced manual wheelchair designed to provide reliable support, comfort, and adjustability for users requiring a higher weight capacity. Engineered with durability in mind, the RubiX<sup>2</sup> XL accommodates users up to 170 kg, making it ideal for everyday mobility in both home and care settings.

With adjustable seat depth, backrest angle, and footplate height, this chair delivers personalised comfort and promotes proper posture. Its heavy-duty frame ensures strength without compromising manoeuvrability, allowing users to move with confidence and independence.

### **Key Features:**

- Heavy-duty frame with a 170 kg weight capacity
- Adjustable seat depth, back angle, and footrest height for tailored comfort
- Reinforced cross-brace and durable construction for long-term reliability
- Swing-away, height-adjustable leg rests and removable armrests
- Large rear wheels and ergonomic hand rims for easy self-propulsion

- Compatible with a range of seating and positioning accessories
- Ideal for home use, rehab centres, and long-term care facilities

## Specifications:

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight
074000-001-52	Breezy Rubix 2 XL - 52cm Wide	410 - 510 mm	52cm (20.4")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg
074000-001-56	Breezy Rubix 2 XL - 56cm Wide	410 - 510 mm	56cm (22")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg
074000-001-60	Breezy Rubix 2 XL - 62cm Wide	410 - 510 mm	62cm (24.4")	410 - 460 mm	375 - 510 mm	170 kg	16.8 kg