



## BRIBIE BARIATRIC FOUR POINT CANE



**SKU:** RBE806

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Bribie Bariatric Four Point Cane** provides enhanced stability and support for users requiring additional balance assistance. Designed with a reinforced steel frame and four-point base, it offers a secure and durable mobility aid suitable for users up to 270 kg.

### Features

- Four-point base design for increased stability and balance
- Adjustable handle height with large double push-buttons for easy adjustment
- Constructed from tough reinforced steel for strength and durability
- Supports users up to 270 kg
- Lightweight design for ease of handling

### Specifications

- **Overall Height:** 71 – 93 cm (adjustable)
- **Base Length:** 24 cm

- **Base Width:** 18 cm
- **Base Height:** 8 cm
- **Unit Weight:** 1.5 kg
- **Weight Capacity:** 270 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)