



BRIBIE BARIATRIC FOUR POINT CANE

\$150.00



SKU: RBE806

PRODUCT DESCRIPTION

Purpose & Benefits

The **Bribie Bariatric Four Point Cane** provides enhanced stability and support for users requiring additional balance assistance. Designed with a reinforced steel frame and four-point base, it offers a secure and durable mobility aid suitable for users up to 270 kg.

Features

- Four-point base design for increased stability and balance
- Adjustable handle height with large double push-buttons for easy adjustment
- Constructed from tough reinforced steel for strength and durability
- Supports users up to 270 kg
- Lightweight design for ease of handling

Specifications

• Overall Height: 71 - 93 cm (adjustable)

• Base Length: 24 cm



• Base Width: 18 cm

• Base Height: 8 cm

• Unit Weight: 1.5 kg

• Weight Capacity: 270 kg