



## **BUCKINGHAM EASY WIPE**

**\$119.00**



**SKU:** PAT-AA2651

## **PRODUCT DESCRIPTION**

## **PURPOSE AND BENEFITS**

The Buckingham Easywipe assists with personal hygiene for those who find reaching difficult and removes the need to touch the paper after using the toilet. A long handled aid that holds toilet paper or wet wipes. It has a lightweight curved plastic shaft. The toilet paper or wipe is attached to a smooth silicon end. A bean shaped recess holds the paper in place for cleaning. A push button at the opposite end releases the toilet paper or wipe after use.

## **FEATURES**

- Easy to use mechanism grips the toilet tissue firmly and then allows it to be released
- A smooth rounded and ergonomic shape enables use for those with limited manual dexterity
- Suitable for use with either toilet tissue or wet wipes
- Supplied with handy carry bag for hygienic storage
- Latex free