




BYACRE CARBON BACKREST



SKU: 20105-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **byACRE Carbon Backrest** provides comfortable lower back support while seated on your rollator, enhancing posture and stability. Designed for easy attachment, it offers a secure fit for compatible byACRE rollator models.

FEATURES

- Provides excellent lower back support while seated on your rollator
- Adjustable height for a personalised fit (up to 25 mm below the rollator's tallest handle position)
- Simple and secure attachment using the built-in Velcro strap and metal bracket system
- Durable, lightweight PE construction for lasting comfort and reliability

This accessory is designed for use with the [byACRE Carbon Ultralight](#) and [byACRE Carbon Overland](#) rollators.

SPECIFICATIONS

Size	Length	Width	Height	Material	Weight
Compact / Regular	85 cm	0.35 cm	15 cm	PE	0.2 kg
Wide Track	89 cm	0.35 cm	15 cm	PE	0.2 kg

Note: Accessory only – walker sold separately.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au