



BYACRE CARBON OVERLAND - LIGHTWEIGHT OFF-ROAD WALKER



SKU: 10038-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF BYACRE CARBON OVERLAND - LIGHTWEIGHT OFF-ROAD WALKER

The byACRE Carbon Overland is the world's lightest all-terrain rollator, weighing in at just 6.7 kg. It is made from strong and durable carbon fibre, making it capable of handling even the most challenging terrain. The Carbon Overland also features ultra-durable pneumatic rubber tyres with maximum shock absorption, rigid front forks and sealed bearings, and a comfortable built-in seat.

- **Freedom of Movement:** Whether you're in a park or on a rough trail, the Carbon Overland ensures you move without limitations.
- **Lightweight Design:** Its carbon fibre build means it's easy to lift, carry, or store without straining oneself.
- **Stylish and Modern:** The sleek design ensures users can confidently use their walker without feeling self-conscious.
- **Durability:** Built to withstand rough terrains and daily wear and tear, ensuring longevity.
- **Safety:** The braking system combined with the sturdy frame provides users with a sense of security.

FEATURES

- **Advanced Carbon Fiber Frame:** The use of carbon fibre ensures the walker is not only lightweight for easy transportation and manoeuvrability but also incredibly strong and durable for off-road conditions.
- **All-Terrain Wheels:** Equipped with large, durable wheels designed to tackle a variety of outdoor terrains, ensuring smooth passage over obstacles and uneven surfaces.

- **Ergonomic Handgrips:** The handgrips are designed for comfort over long periods of use, reducing the risk of fatigue and providing a secure grip in all weather conditions.
- **Easy-Fold Mechanism:** Allows for quick and easy folding of the walker for storage or transport, making it convenient for adventures far from home.
- **Adjustable Height Settings:** Ensures a personalised fit for users of different heights, improving posture and comfort during use.

ATTACHMENTS

The byACRE Carbon Overland has a range of attachments available.

[Click here](#) to view all of our walker attachments.

SPECIFICATIONS

- **Lifetime warranty on carbon fibre frame** | 2-year warranty on all other parts*
- **Material:** High-grade carbon fibre
- **Colour Options:** Overland Brown, Overland Green

Weight	6.7 kg
Handle height	82-95 cm
Length	73 cm
Seat width	47 cm
Seat height	61 cm
Width	66 cm
Width when folded	25 cm
Max. user weight	120kg
Wheel diameter	25 cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au