



## BYACRE CARBON OVERLAND - LIGHTWEIGHT OFF-ROAD WALKER

**\$1,194.00** Sale Available Online Only!



**SKU:** 10038-H

### PRODUCT DESCRIPTION

#### PURPOSE & BENEFITS OF BYACRE CARBON OVERLAND - LIGHTWEIGHT OFF-ROAD WALKER

The byACRE Carbon Overland is the world's lightest all-terrain rollator, weighing in at just 6.7 kg. It is made from strong and durable carbon fibre, making it capable of handling even the most challenging terrain. The Carbon Overland also features ultra-durable pneumatic rubber tyres with maximum shock absorption, rigid front forks and sealed bearings, and a comfortable built-in seat.

- **Freedom of Movement:** Whether you're in a park or on a rough trail, the Carbon Overland ensures you move without limitations.
- **Lightweight Design:** Its carbon fibre build means it's easy to lift, carry, or store without straining oneself.
- **Stylish and Modern:** The sleek design ensures users can confidently use their walker without feeling self-conscious.
- **Durability:** Built to withstand rough terrains and daily wear and tear, ensuring longevity.
- **Safety:** The braking system combined with the sturdy frame provides users with a sense of security.

### FEATURES

- **Advanced Carbon Fiber Frame:** The use of carbon fibre ensures the walker is not only lightweight for easy transportation and maneuverability but also incredibly strong and durable for off-road conditions.
- **All-Terrain Wheels:** Equipped with large, durable wheels designed to tackle a variety of outdoor terrains, ensuring smooth passage over obstacles and uneven surfaces.

- **Ergonomic Handgrips:** The handgrips are designed for comfort over long periods of use, reducing the risk of fatigue and providing a secure grip in all weather conditions.
- **Easy-Fold Mechanism:** Allows for quick and easy folding of the walker for storage or transport, making it convenient for adventures far from home.
- **Adjustable Height Settings:** Ensures a personalized fit for users of different heights, improving posture and comfort during use.

## ATTACHMENTS

The byACRE Carbon Overland has a range of attachments available.

[Click here](#) to view all of our walker attachments.

## SPECIFICATIONS

- **Lifetime warranty on carbon fibre frame** | 2-year warranty on all other parts\*
- **Material:** High-grade carbon fibre
- **Color Options:** Overland Brown, Overland Green

|                          |          |
|--------------------------|----------|
| <b>Weight</b>            | 6.7 kg   |
| <b>Handle height</b>     | 82-95 cm |
| <b>Length</b>            | 73 cm    |
| <b>Seat width</b>        | 47 cm    |
| <b>Seat height</b>       | 61 cm    |
| <b>Width</b>             | 66 cm    |
| <b>Width when folded</b> | 25 cm    |
| <b>Max. user weight</b>  | 120kg    |
| <b>Wheel diameter</b>    | 25 cm    |