

## BYACRE CARBON ULTRALIGHT WALKER

**\$1,194.00** Sale Available Online Only!



**SKU:** 2110-H



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS OF THE BYACRE CARBON ULTRALIGHT WALKER

The byACRE Carbon Ultralight Walker is engineered to enhance mobility and independence for individuals requiring support while walking. This revolutionary product stands out for its exceptional blend of strength, durability, and lightweight design, making it an ideal companion for daily use, travel, and outdoor activities. It offers an unmatched level of comfort and ease of use, enabling users to navigate various terrains with confidence and ease. The walker's sleek and modern design also ensures users can enjoy functionality without compromising on style.

### FEATURES

- **Ultra-lightweight Design:** Crafted from advanced carbon fibre material, this walker is both incredibly strong and surprisingly light, making it easy to carry, maneuver, and store.
- **Ergonomic Handles:** The ergonomically designed handles provide a comfortable grip, reducing hand fatigue and enhancing user comfort during extended periods of use.
- **Adjustable Height:** With simple, tool-free adjustments, the walker can be tailored to the perfect height, ensuring an optimal and safe walking posture for users of different statures.
- **Easy-to-Use Braking System:** Features an intuitive braking system that is easily accessible, offering reliable stopping power and control for added safety.
- **Foldable Design for Easy Storage and Transport:** The walker folds effortlessly in seconds, making it convenient for storage in small spaces and transport in a car or on public transportation.
- **Durable, All-Terrain Wheels:** Equipped with robust wheels designed to handle various surfaces, from indoor floors to outdoor pavements, ensuring smooth mobility across different environments.

## SPECIFICATIONS

- **Material:** High-quality carbon fibre
- **Safe Working Load:** 130kg
- **Colour Options:** Carbon Black, Oyster White, Strawberry Red

| Technical Data        | Compact | Regular | Wide Track |
|-----------------------|---------|---------|------------|
| Weight                | 4.8kg   | 4.8kg   | 4.95kg     |
| Handle Height         | 79-91cm | 79-91cm | 79-91cm    |
| Width between handles | 38cm    | 39cm    | 44cm       |
| Seat width            | 41cm    | 42cm    | 47cm       |
| Seat height           | 56cm    | 61cm    | 61cm       |
| Width                 | 61.5cm  | 62.5cm  | 66.5cm     |
| Width when folded     | 25cm    | 25cm    | 25cm       |
| Length                | 66cm    | 66cm    | 66cm       |

## ATTACHMENTS

The byACRE Ultralight has a range of attachments available.

[Click here](#) to view all of our walker attachments.

## VIDEO

<https://youtu.be/XBsf2j0q0c8>