

BYACRE CARBON ULTRALIGHT WALKER



SKU: 2110-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF THE BYACRE CARBON ULTRALIGHT WALKER

The byACRE Carbon Ultralight Walker is engineered to enhance mobility and independence for individuals requiring support while walking. This revolutionary product stands out for its exceptional blend of strength, durability, and lightweight design, making it an ideal companion for daily use, travel, and outdoor activities. It offers an unmatched level of comfort and ease of use, enabling users to navigate various terrains with confidence and ease. The walker's sleek and modern design also ensures users can enjoy functionality without compromising on style.

FEATURES

- **Ultra-lightweight Design:** Crafted from advanced carbon fibre material, this walker is both incredibly strong and surprisingly light, making it easy to carry, manoeuvre, and store.
- **Ergonomic Handles:** The ergonomically designed handles provide a comfortable grip, reducing hand fatigue and enhancing user comfort during extended periods of use.
- **Adjustable Height:** With simple, tool-free adjustments, the walker can be tailored to the perfect height, ensuring an optimal and safe walking posture for users of different statures.
- **Easy-to-Use Braking System:** Features an intuitive braking system that is easily accessible, offering reliable stopping power and control for added safety.
- **Foldable Design for Easy Storage and Transport:** The walker folds effortlessly in seconds, making it convenient for storage in small spaces and transport in a car or on public transportation.
- **Durable, All-Terrain Wheels:** Equipped with robust wheels designed to handle various surfaces, from indoor floors to outdoor pavements, ensuring smooth mobility across different environments.

SPECIFICATIONS

- **Material:** High-quality carbon fibre
- **Safe Working Load:** 130kg
- **Colour Options:** Carbon Black, Oyster White, Strawberry Red

Technical Data	Compact	Regular	Wide Track
Weight	4.8kg	4.8kg	4.95kg
Handle Height	79-91cm	79-91cm	79-91cm
Width between handles	38cm	39cm	44cm
Seat width	41cm	42cm	47cm
Seat height	56cm	61cm	61cm
Width	61.5cm	62.5cm	66.5cm
Width when folded	25cm	25cm	25cm
Length	66cm	66cm	66cm

ATTACHMENTS

The byACRE Ultralight has a range of attachments available.

[Click here](#) to view all of our walker attachments.

VIDEO

<https://youtu.be/XBsf2j0q0c8>

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au