



## BYACRE NORDIC PIONEER ROLLATOR



**SKU:** 24102

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **byACRE Nordic Pioneer** is a next-generation rollator designed to combine innovation, sustainability, and superior functionality. Winner of the 2024 German Design Award, it redefines mobility with its lightweight 5.9 kg EPP (Expanded Polypropylene) frame reinforced with an aluminium core. Built for both indoor and outdoor use, it offers smooth handling, unmatched portability, and a premium walking experience that promotes upright posture and independence.

### FEATURES

- Winner of the **German Design Award 2024** for innovation and usability
- Lightweight at only **5.9 kg** for effortless portability and lifting
- Constructed from eco-friendly, durable EPP foam with an **aluminium core** for strength and longevity
- **Soft “bump-proof” frame** provides safety and shock absorption on uneven surfaces
- **Large, low-friction wheels** for smooth rolling on mixed terrains
- **Forward-facing ergonomic handles** promote upright posture and reduce strain on shoulders and back
- **Easy folding mechanism** folds and unfolds in seconds for convenient transport and storage
- **Built-in seat** with stain-resistant and water-repellent material for rest anywhere, anytime
- **Adjustable handle height** (84–94 cm) with simple push-button adjustment for personalised comfort
- **High-quality PUR tyres** offer excellent traction, shock absorption, and are puncture-proof
- Sealed glide bearings ensure quiet, wobble-free performance
- Stylish **Glacier Green** finish for a modern, sophisticated look

- Delivered fully assembled and ready to use

## SPECIFICATIONS

- **Brand:** byACRE
- **Model:** Nordic Pioneer
- **Award:** German Design Award Winner 2024
- **Weight:** 5.9 kg
- **Handle Height:** 84 – 94 cm (4 height settings)
- **Width Between Handles:** 43 – 43.5 cm
- **Grip to Brake Distance:** 6.2 cm
- **Seat Height:** 61.5 cm
- **Seat Width:** 47 cm
- **Maximum User Weight:** 100 kg
- **Wheel Diameter:** 22.5 cm
- **Width (Folded):** 24 cm
- **Length:** 69 cm
- **Material:** Eco-friendly EPP outer with aluminium core
- **Colour:** Glacier Green

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)