



BYACRE OVERLAND BACK REST

\$99.00 Original price was: \$99.00.<u>\$89.10</u>Current price

is: \$89.10.



SKU: 20206

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **byACRE Overland Back Rest** provides comfortable lower-back support while seated on your rollator. Designed to replicate the feel of a chair back, it offers both comfort and stability during rest breaks, making it an ideal companion for longer strolls or outdoor adventures.

FEATURES

- Supports the lower back while seated on your rollator
- Ergonomic shape designed to simulate a traditional chair back
- Height adjustable for a personalised and comfortable fit
- Made from durable, lightweight PE thermoplastic
- Soft, stable and weather-resistant construction
- Easy to attach and detach without tools
- Compatible with the byACRE Overland and Wide Track byACRE Ultralight

SPECIFICATIONS

Specification	Detail
Compatibility	byACRE Overland and Wide Track byACRE Ultralight
Material	PE (lightweight, durable thermoplastic)
Colour	Dark grey
Weight	0.2 kg
Dimensions	Width: 89 cm Height: 15 cm



Note: Accessory part only – rollator not included.