




## BYACRE SCANDINAVIAN BUTLER INDOOR WALKER

~~\$995.00~~ Original price was: \$995.00. ~~\$922.50~~ Current price is: \$922.50. Sale Available Online Only!



SKU: 23101

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The byACRE Scandinavian Butler Indoor Walker is an elegantly designed, highly functional mobility aid that combines the practicality of a traditional walker with the aesthetic appeal of Scandinavian design. Created to enhance indoor mobility, this walker is perfect for individuals who require assistance with walking but also value style and design in every aspect of their lives.

### FEATURES

- **Sleek Scandinavian Design:** The byACRE Scandinavian Butler Indoor Walker features a minimalist design that is both stylish and functional. Its clean lines and modern aesthetics blend seamlessly with any home decor, making it a discrete and attractive mobility aid.
- **Lightweight Construction:** Crafted from high-quality, lightweight materials, this walker is easy to maneuver around the home, enabling users to maintain their independence without the bulkiness of traditional walkers.
- **Built-In Tray and Basket:** Equipped with a detachable tray and a basket, the Scandinavian Butler allows users to easily transport items such as drinks, meals, and personal belongings from one room to another, facilitating independence and convenience.
- **Adjustable Handles:** The handles are adjustable, ensuring a comfortable grip and proper posture while using the walker. This feature allows for customization to fit the user's height and preferences.
- **Sturdy and Stable:** Despite its lightweight design, the walker provides robust support and stability, thanks to its well-engineered frame and durable wheels, ensuring users feel secure at all times.
- **Brake System:** Features an easy-to-use brake system that enhances safety by allowing the

user to control the speed and secure the walker in place when needed.

## **SPECIFICATIONS**

- Colour **Oyster White**
- Weight **10kg**
- Max user weight **150kg**
- Length **63.3cm**
- Height **83.5cm**
- Width **58.8cm**
- Handle height **83.5-97.5cm**
- Handle width **43cm**
- Boxed dimensions: **57.6cm x 49.7cm x 72cm - 12kg**
- Material: **Steel frame, aluminium, POM and PA**
- Wheels: **ABS coated with soft rubber**
- Wheels diameter: **18cm**