



CAMEL LIFTING CUSHION WITH AIRFLO 24 | SIT TO STAND AID

\$9,392.00



SKU: LSS403810

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Camel Lifting Cushion is an inflatable device that aids in lifting individuals who have fallen. It can be used indoors or outdoors, providing a dignified lift while ensuring the safety of both the caregiver and the person being lifted. This innovative product is especially suitable for larger or bariatric individuals, with a maximum weight capacity of 320 kg. It comes equipped with the Airflo 24 compressor, which operates using rechargeable batteries.

The standout feature of this cushion is its ability to allow a single person to safely raise a fallen client to a secure position, although it is recommended to have two caregivers present. We offer comprehensive training and demonstrations to ensure proper usage. Feel free to contact us for more information.

Who can benefit from the Camel Lifting Cushion? The Camel Lifting Cushion is an inflatable device specifically designed to assist in lifting fallen individuals. It is capable of supporting individuals weighing up to 320 kg, effortlessly raising them from the floor to a seated position. Additionally, it provides excellent back support for individuals who require additional stability during the lifting process. By utilizing the inflatable Camel cushion and Airflo 24 air compressor, the physical strain on caregivers associated with manual handling is significantly reduced.

FEATURES

- The Camel Lifting Cushion is particularly suitable for very large or bariatric people (lifts up to 320kg)
- Lightweight
- Simple to use
- Appropriate for any age

- Portable
- Packs away into a compact and discreet bag
- For use indoors or out
- Minimises moving & handling risks
- Enhances independent living
- Battery operated
- Airflo controlled

SPECIFICATIONS

- Max user weight: 320kgs
- Inflated Seat Height: 21.5"
- Inflated Width: 28"
- Inflated backrest width: 24.5"
- Inflated seat depth: 19"
- Inflated backrest height: 32.5"
- Overall inflated length: 44"
- Deflated rolled dimensions: 30.5" x 8"
- Weight: 6.3kg

[Click to view the Camel Lifting Cushion with Airflo 24 Brochure](#)