



ALRICK 3250 SERIES CHANGE TABLE / EXAMINATION COUCH



SKU: 3250

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Alrick 3250 Series Changetable/Examination Couch Series has taken the place of both the 3006 and the 3008 Series Changetable. Its new features make it the ideal option for multiple environments with a 450 to 900 mm height adjustment range and a massive 250 kg safe working load. The as well as the lift function is a vertical lift and for added flexibility the backrest is still available as an option on this changetable (3250/2). Central locking castors are also now an option on the 3250 Changetable.

FEATURES

- 3250 Series Changetable / Examination Couch
- Electric Hi Low
- Handset
- Standard 125mm Braking Castors
- Charcoal Padded Vinyl Top
- Two Tone Powder Coated Frame (APO Grey / Ironstone) Overall Dimensions 1880mm x 700mm
- Height Range 450 - 900mm
- Weight Capacity 250kg
- Replaces 3006 & 3008 Series

SPECIFICATIONS

Minimum Height - 450mm

Maximum Height - 900mm

Overall Width – 700mm
Overall Length – 1880mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au