




## **ALRICK 3250 SERIES CHANGE TABLE / EXAMINATION COUCH**



**SKU: 3250**

### **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Alrick 3250 Series Changetable/Examination Couch Series has taken the place of both the 3006 and the 3008 Series Changetable. Its new features make it the ideal option for multiple environments with a 450 to 900 mm height adjustment range and a massive 250 kg safe working load. The as well as the lift function is a vertical lift and for added flexibility the backrest is still available as an option on this changetable (3250/2). Central locking castors are also now an option on the 3250 Changetable.

### **FEATURES**

- 3250 Series Changetable / Examination Couch
- Electric Hi Low
- Handset
- Standard 125mm Braking Castors
- Charcoal Padded Vinyl Top
- Two Tone Powder Coated Frame (APO Grey / Ironstone) Overall Dimensions 1880mm x 700mm
- Height Range 450 - 900mm
- Weight Capacity 250kg
- Replaces 3006 & 3008 Series

### **SPECIFICATIONS**

Minimum Height - 450mm

Maximum Height - 900mm

Overall Width – 700mm

Overall Length – 1880mm