



## **CAREQUIP HEALTHWEIGH CHAIR SCALES**



**SKU:** JE0660

### **PRODUCT DESCRIPTION**

## Purpose & Benefits of Carequip Healthweigh Chair Scales

The Carequip Healthweigh Chair Scales are expertly crafted to provide precise and reliable weight measurements for individuals who require seated weighing options, particularly in medical and care settings. These scales are designed to accommodate patients comfortably, ensuring their safety and stability during the weighing process. They are ideal for use in hospitals, clinics, and elderly care facilities, offering an essential tool for monitoring patient health.

#### **Features**

- **Ergonomic Design**: The chair is contoured to fit the shape of the patient, enhancing comfort during use.
- Adjustable Components: Features folding armrests and footrests that can be raised independently, adapting to the needs of each patient for easier access and positioning.
- **Movement Compensation Technology**: Helps achieve accurate weight measurements even if the patient is slightly moving, ensuring reliability.
- **BMI Functionality**: Includes a Body Mass Index calculation feature, adding valuable health metric assessments for patient care management.
- **User-Friendly Interface**: Equipped with a 5-digit LCD screen that provides clear and immediate weight readings.
- **Robust Construction**: Built with durable materials to withstand frequent use in professional healthcare environments.
- **Hygienic Surface**: The surfaces are easy to clean, maintaining hygiene standards required in medical settings.



# **Specifications**

• Display: 5 Digit LCD screen

• **Graduations**: Measures in increments of 0.1 kg for precise weight tracking.

• Functions: Tare/Zero/Hold/kg/lb options

• Safe Working Load: 300 kg

Dimensions:

Seat Width: 430 mm
Seat Depth: 420 mm
Overall Width: 630 mm
Overall Depth: 850 mm
Overall Height: 965 mm

Weight: 18 kg

The Carequip Healthweigh Chair Scales combine functionality, comfort, and precision, making them a critical component in the effective management of patient care. Their robust design and advanced features ensure that they meet the demands of today's healthcare industry, prioritizing patient care, safety, and comfort.

For further information go to medimart.com.au or give us a call: 1300 633 063

#### SA

#### **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

#### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

## Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

#### **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

#### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



# NT

## **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 <a href="mailto:support@medimarnt.com.au">support@medimarnt.com.au</a>