



## CAREQUIP MOBILE IV POLE

**\$210.00**



**SKU: GB0010**

## PRODUCT DESCRIPTION

### Purpose & Benefits

The CareQuip Mobile IV Pole provides a practical and dependable solution for medical settings where fluid administration needs to be efficient and mobile. Designed with user convenience in mind, this pole features a stable base and smooth-rolling castors that allow effortless movement between rooms or patient bedsides.

With adjustable height and twin hooks, it accommodates a range of IV bags securely and is ideal for use in clinics, hospitals, or at home. Its durable construction ensures long-lasting performance while offering the flexibility to support diverse patient care needs.

Keep patients comfortable and monitored without needing to move large equipment

Easily relocate the pole thanks to its four smooth-glide castors

Designed for stable performance, even during frequent repositioning

Ideal for home care, aged care, and clinical use

### Features

- Mobile IV Pole with Twin Bag Hooks
- Four 50 mm Castors for Easy Mobility
- 5-Point Star Base for Enhanced Stability

- Durable Powder-Coated Steel Construction
- Lightweight & Easy to Handle
- Designed for Hospitals, Aged Care & Home Use
- Australian Standard Quality

## **Specifications**

- Overall Height: 1240-2130 mm
- Base Diameter: 560 mm