

## CARICHAIR PUMP WITH MOBI PLUS CUSHION

**\$1,960.00 - \$2,053.00** Price range: \$1,960.00 through \$2,053.00



**SKU:** APC-CARI-MO-PL-H



## PRODUCT DESCRIPTION

### Purpose & Benefits of the CariChair Pump with Mobi Plus Cushion

The CariChair Pump with Mobi Plus Cushion is designed to provide effective pressure relief and comfort for individuals at risk of developing pressure sores while seated. This system combines an alternating pressure pump with a cushion to promote healthy blood circulation and reduce prolonged pressure on specific areas, enhancing user comfort and minimising the risk of pressure ulcers.

### Features

- **Alternating Pressure Technology** – The pump system alternates pressure within the cushion to alleviate and redistribute pressure points, reducing the risk of pressure sores.
- **Compact and Lightweight Design** – The cushion measures 450 mm in length, 360 mm in width, and 90 mm in height, making it suitable for various seating arrangements without adding significant bulk.
- **User-Friendly Operation** – The system is designed for easy setup and operation, ensuring that users or caregivers can manage it without difficulty.
- **Quiet Pump Mechanism** – The pump operates quietly, ensuring minimal disruption during use.

### Specifications

- **Dimensions:** 450 mm (L) x 360 mm (W) x 90 mm (H)
- **Weight:** 1.5 kg
- **Maximum User Weight Capacity:** 150 kg
- **Risk Category:** Suitable for individuals at medium to high risk of pressure injury

The CariChair Pump with Mobi Plus Cushion offers an effective solution for individuals seeking to

prevent pressure sores and enhance comfort during extended periods of sitting.