


CARICHAIR PUMP WITH MOBI PLUS CUSHION



SKU: APC-CARI-MO-PL-H



PRODUCT DESCRIPTION

Purpose & Benefits of the CariChair Pump with Mobi Plus Cushion

The CariChair Pump with Mobi Plus Cushion is designed to provide effective pressure relief and comfort for individuals at risk of developing pressure sores while seated. This system combines an alternating pressure pump with a cushion to promote healthy blood circulation and reduce prolonged pressure on specific areas, enhancing user comfort and minimising the risk of pressure ulcers.

Features

- **Alternating Pressure Technology** – The pump system alternates pressure within the cushion to alleviate and redistribute pressure points, reducing the risk of pressure sores.
- **Compact and Lightweight Design** – The cushion measures 450 mm in length, 360 mm in width, and 90 mm in height, making it suitable for various seating arrangements without adding significant bulk.
- **User-Friendly Operation** – The system is designed for easy setup and operation, ensuring that users or caregivers can manage it without difficulty.
- **Quiet Pump Mechanism** – The pump operates quietly, ensuring minimal disruption during use.

Specifications

- **Dimensions:** 450 mm (L) x 360 mm (W) x 90 mm (H)
- **Weight:** 1.5 kg
- **Maximum User Weight Capacity:** 150 kg
- **Risk Category:** Suitable for individuals at medium to high risk of pressure injury

The CariChair Pump with Mobi Plus Cushion offers an effective solution for individuals seeking to

prevent pressure sores and enhance comfort during extended periods of sitting.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au