



CARILEX COZINY 300 MATTRESS



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Carilex Coziny 300 Mattress

The **Carilex Coziny 300 Mattress** is a premium static foam support surface designed to deliver superior comfort and pressure redistribution for individuals at medium to high risk of pressure injuries. With its tri-layer foam construction and contoured, castellated top layer, the Coziny 300 promotes effective immersion, envelopment, and airflow—minimising shear and friction on vulnerable skin areas.

Ideal for use in hospitals, aged care facilities, and community settings, this mattress offers enhanced postural support and durability. Its vapour-permeable, waterproof cover ensures optimal hygiene while supporting skin integrity, making it a reliable solution for long-term patient care.

Key Features

- Tri-layer foam construction for advanced pressure redistribution
- Castellated top layer enhances immersion, airflow, and comfort
- Suitable for patients at medium to high risk of pressure injuries
- Vapour-permeable, waterproof, and easy-clean cover
- Promotes spinal alignment and reduces shear forces

- Durable, lightweight design for ease of handling and transport
- Ideal for hospital, aged care, and home care environments

Specifications

Product Name	Coziny 300 Cot
Product Category	Mattress Replacement
Risk	Med-High
Waterlow Scale	15+
Pressure Sore Grading	3
Safe Working Load	25kg
Dynamic Mode	✓
Static Mode - Constant Low Pressure	✓
Intelligent Pressure Sensor	✓
Low Pressure/Power Failure Alarms + Battery Backup	✓
Number of Air Cells	39 x 1.5"
CPR Fast Deflation	via pump
Cover - 2 Way Stretch	-
Dartex Cover - 4 Way Stretch	✓
Power Cable Management System	-
Washable Cover	✓
Digital Control	✓
True Low Air Loss	-
Integrated Foam or Air Base	Foam
Cycle Times (mins)	10, 15, 20
Upright Mode	-
Firm Mode	✓
Memory Recall	✓
Auto Panel Lock	✓
Battery Backup for Alternating Mode	12 hours
Mattress Size - Single	-
Mattress Size - Euro King Single	-
Mattress Size - King Single	-
Mattress Size - Extra Wide	-
Mattress Size - Double	-
Mattress Size - Queen	-

Mattress Size - Other

70 x 140 x 8cm

Mattress Firm Side Bolsters

–