

CARILEX DUALPLUS APM MATTRESS



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Dual Plus mattress is customisable to suit individual patient needs. Localised healing can be created by easily swapping any of the 8" cells with the 5" bridging cell, providing reduced pressure zones, for example, heel zone.

SPECIFICATIONS

Bed Size	Single
Max User Weight	250kg
Overall Height	210mm
Overall Length	2000mm
Overall Width	820mm
Product Weight	11.3kg