


## CARILEX DUALPLUS APM MATTRESS



SKU: N/A

### PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Dual Plus mattress is customisable to suit individual patient needs. Localised healing can be created by easily swapping any of the 8" cells with the 5" bridging cell, providing reduced pressure zones, for example, heel zone.

### SPECIFICATIONS

<b>Bed Size</b>	Single
<b>Max User Weight</b>	250kg
<b>Overall Height</b>	210mm
<b>Overall Length</b>	2000mm
<b>Overall Width</b>	820mm
<b>Product Weight</b>	11.3kg