

CARILEX DUALPLUS PUMP



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Carilex DualPlus Pump

The **Carilex DualPlus Pump** is a compact, high-efficiency digital control unit designed to operate Carilex dynamic mattress systems with precision and reliability. Supporting both alternating pressure and continuous low pressure modes, the DualPlus Pump delivers clinically effective pressure relief for patients at high to very high risk of pressure injuries.

Engineered with user-friendly digital controls and automatic weight-sensing technology, the pump adapts in real time to ensure consistent pressure redistribution. Its quiet operation and integrated safety features make it ideal for use in hospitals, aged care facilities, and home care settings where both performance and comfort are critical.

Key Features

- Supports alternating pressure and continuous low pressure therapy
- Compatible with a range of Carilex dynamic mattress systems
- Auto weight-sensing for responsive, individualised pressure adjustment
- Digital interface with simple, intuitive controls
- Quiet operation for patient comfort in all care environments

- Integrated safety features including visual and audible alarms
- Designed for high to very high-risk pressure injury prevention

Specifications

TECHNICAL DETAILS

Risk/Category	High to Very High
Waterlow Score	20+
Pressure Sore Grading	Up to 4
Safe Working Load	30 to 250kg
Compatible Support Surfaces	Dual Plus Modular Mattress – Single and King Single
Travel / Weight	Pump weighs 2.3kg
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Warranty	2 years