

CARILEX ENTRIX 7 SINGLE MATTRESS



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Carilex Entrix 7 Single Mattress

The **Carilex Entrix 7 Single Mattress** is a dynamic pressure care mattress replacement system designed to provide advanced support and pressure redistribution for individuals at high to very high risk of pressure injuries. Combining alternating pressure and continuous low pressure therapy modes, the Entrix 7 ensures optimal skin protection and comfort across a variety of care settings.

Ideal for use in hospitals, aged care facilities, and home environments, the system features an intuitive digital pump, low-noise operation, and auto-weight sensing to deliver precise pressure adjustments. The multi-layered cell-on-cell design provides stability and protection during power failure, while the breathable, waterproof cover enhances hygiene and ease of maintenance.

Key Features

- Alternating and continuous low pressure modes for individualised therapy
- Cell-on-cell air mattress design with integrated static base for safety during power failure
- Digital pump with auto weight-sensing for accurate pressure adjustment
- Low noise and vibration for patient comfort
- Vapour-permeable, waterproof, and easy-clean cover

- Suitable for high to very high-risk patients
- Ideal for hospital, aged care, and in-home care settings

Specifications

- Size: Single 90 x 200 x 18 cm; King Single 105 x 200 x 18 cm;
- Safe Working Load: 150kg