

## CARILEX ENTRIX NX PUMP WITH 5" OVERLAY



SKU: N/A

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

- 5" mattress overlay system designed to sit on top of your existing mattress.
- Envelops and off loads patient body mass by expanding the therapy contact area and depth.
- Easy connection between power unit and mattress.
- No fuss cable and cord management.
- CPR pull tab deflates the entire mattress in seconds.
- 2 way stretch washable cover.
- Available in single and king single.

### SPECIFICATIONS

Overall Length	2000 mm
Overall Width	900 mm
Thickness	130 mm
Pressure Sore Grading	2
Waterlow Scale	15+
Product Weight	4.1 kg
Max User Weight	150 kg
Brand	Carilex