


CARILEX LUX MATTRESS



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Carilex Lux Mattress

The **Carilex Lux Mattress** is a high-performance dynamic mattress replacement system designed to deliver optimal pressure redistribution and advanced patient support. Intended for individuals at high to very high risk of pressure injuries, the Lux Mattress combines intelligent therapy modes with premium comfort features to meet the demands of acute, aged care, and home environments.

Featuring alternating pressure, continuous low pressure, and max firm modes, the Lux system adapts to individual needs while promoting circulation and reducing pressure-related complications. Its multi-layered cell-on-cell design ensures enhanced safety during power failure, while the low-noise digital pump and vapour-permeable cover contribute to patient comfort and infection control.

Key Features

- Dynamic mattress replacement system for high to very high-risk patients
- Multiple therapy modes, including alternating pressure and CLP
- Cell-on-cell construction for support during power failure
- Digital pump with user-friendly controls and auto-weight adjustment
- Low-noise operation ideal for 24-hour care environments

- Vapour-permeable, waterproof cover for hygiene and skin protection
- Suitable for hospitals, aged care, and in-home care settings

Specifications

TECHNICAL DETAILS

Risk/Category	Very High
Waterlow Score	20+
Pressure Sore Grading	4
Safe Working Load	Up to 180kg
Air Cell Configuration	20 x 7" cells
Cycle Times	10 minutes, 2:1
Available Size / Dimensions	Single 90cm x 200cm x 18cm King Single 107cm x 200cm x 18cm
Securing Mattress to Bed	Buckle straps
Travel / Weight	Mattress rolls up into a bag, weighing 7kg
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Compatible Pumps	Entrix NX and Centrius
Warranty	2 years

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au