

## **CARILEX THERAFLO AP MATTRESS**





SKU: N/A

## **PRODUCT DESCRIPTION**

## **Purpose & Benefits of Carilex TheraFlo AP Mattress**

The **Carilex TheraFlo AP Mattress** is a non-powered, air and foam hybrid mattress system designed to provide effective pressure redistribution for patients at medium to high risk of pressure injuries. With its self-adjusting air cell system, the TheraFlo AP delivers dynamic support without the need for an external pump or power source, making it ideal for environments where simplicity, quiet operation, and reliability are essential.

Each air cell automatically responds to the user's weight and movement, providing consistent immersion and envelopment to reduce shear and pressure forces. This low-maintenance solution is suitable for hospitals, aged care facilities, and home care, supporting both patient comfort and carer efficiency.

# **Key Features**

- Non-powered air and foam hybrid mattress system
- Self-regulating air cell technology for automatic pressure adjustment
- No pump or electricity required silent operation
- Consistent pressure redistribution for medium to high-risk patients
- Lightweight and easy to set up and maintain



- Vapour-permeable, waterproof cover for hygiene and comfort
- Suitable for use in hospital, aged care, and home care settings

## **Specifications**

#### **TECHNICAL DETAILS**

Risk/Category Very High

Waterlow Score 20+ Pressure Sore Grading 4

Safe Working Load Up to 280kg Air Cell Configuration 20 x 8" cells

Cycle Times 5, 10, 15 minutes

90cm x 200cm x 25cm 99cm x 200cm x 25cm 107cm x 200cm x 25cm

Available Size / Dimensions

128cm x 200cm x 25cm 137cm x 200cm x 25cm 152cm x 200cm x 25cm

Securing Mattress to Bed Buckle straps

Clinical Settings Suitable for Hospital, Home Care and Nursing Home environments

Compatible Pumps Theraflo AP Warranty 2 years

For further information go to medimart.com.au or give us a call: 1300 633 063

#### SA

Brighton ShowroomElizabeth Vale ShowroomReynella Showroom343 Brighton Rd, NorthPlayford Health Hub, 44 John RiceShop 1/141 Sherriffs Rd,Brighton SA 5048Ave, Elizabeth Vale SA 5112Reynella SA 516108 7079 877308 7079 877708 7079 8771brighton@medimart.com.auelizabethvale@medimart.com.aureynella@medimart.com.au



## **VIC**

## **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

## **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

## NT

## **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au