

CARILEX THERAFLO AP MATTRESS



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Carilex TheraFlo AP Mattress

The **Carilex TheraFlo AP Mattress** is a non-powered, air and foam hybrid mattress system designed to provide effective pressure redistribution for patients at medium to high risk of pressure injuries. With its self-adjusting air cell system, the TheraFlo AP delivers dynamic support without the need for an external pump or power source, making it ideal for environments where simplicity, quiet operation, and reliability are essential.

Each air cell automatically responds to the user's weight and movement, providing consistent immersion and envelopment to reduce shear and pressure forces. This low-maintenance solution is suitable for hospitals, aged care facilities, and home care, supporting both patient comfort and carer efficiency.

Key Features

- Non-powered air and foam hybrid mattress system
- Self-regulating air cell technology for automatic pressure adjustment
- No pump or electricity required – silent operation
- Consistent pressure redistribution for medium to high-risk patients
- Lightweight and easy to set up and maintain

- Vapour-permeable, waterproof cover for hygiene and comfort
- Suitable for use in hospital, aged care, and home care settings

Specifications

TECHNICAL DETAILS

Risk/Category	Very High
Waterlow Score	20+
Pressure Sore Grading	4
Safe Working Load	Up to 280kg
Air Cell Configuration	20 x 8" cells
Cycle Times	5, 10, 15 minutes
Available Size / Dimensions	90cm x 200cm x 25cm
	99cm x 200cm x 25cm
	107cm x 200cm x 25cm
	128cm x 200cm x 25cm
	137cm x 200cm x 25cm
	152cm x 200cm x 25cm
Securing Mattress to Bed	Buckle straps
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Compatible Pumps	Theraflo AP
Warranty	2 years