



## SEAL-TIGHT CAST KNEE PROTECTOR



**SKU:** ST-20324-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of Seal-Tight Cast Knee Protector

The Seal-Tight Cast Knee Protector is specifically designed to safeguard knee casts, bandages, or dressings from water exposure during activities such as bathing, showering, or engaging in water-based activities. This protector helps maintain the dryness and integrity of knee treatments, thereby preventing complications such as infections or delayed healing due to moisture penetration. It is especially useful for individuals recovering from knee surgeries, injuries, or those needing to protect their knee dressings from getting wet.

### Features

- **Durable Waterproof Material:** Constructed from high-grade polyvinyl or polyurethane, providing a strong and reliable waterproof barrier that is resistant to tears and punctures.
- **Easy to Use:** Features a diaphragm seal design that offers a tight fit without the need for additional straps or adhesives, making it easy to put on and remove.
- **Latex-Free:** Suitable for individuals with latex allergies, this protector is made from latex-free materials to prevent skin irritation.
- **Reusable:** Designed for multiple uses, making it a cost-effective solution for patients who require long-term protection for their knee dressings.
- **Comfortable Fit:** The protector is lightweight and designed to offer a snug fit, minimizing discomfort and allowing for ease of movement.
- **Versatile Sizes:** Available in various sizes to accommodate different knee sizes and cast dimensions, ensuring a proper fit for both adults and children.

## Specifications

- **Material:** Durable polyvinyl or polyurethane, latex-free.
- **Size:** Small 41-51cm; Medium 48-58cm; Large 56-66cm;
- **Care Instructions:** The protector can be easily cleaned with a damp cloth and air-dried; it is recommended to store it away from direct sunlight and sharp objects to maintain its integrity.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)