



## **AAXIS 2000ML URIMAX DRAIN BAG A2 - PUSH-PULL DRAIN**

**\$3.00**



**SKU: 12001105**

## **PRODUCT DESCRIPTION**

Urimaax Drainage Bag B/Drain is a latex-free, non-sterile 2000 ml drainage bag with a bottom T-bar drain and a universal stepped catheter connector. Urimaax Drainage Bag B/Drain T-Bar has a graduated front with white frosted backing to allow easy visual readings and hanger holes fit universal hangers.