

BIOCATH CATHETERS - ALL TYPES - BARD

\$15.95 - \$24.00 Price range: \$15.95 through \$24.00



SKU: 10100415



PRODUCT DESCRIPTION

FEATURES

- Flexible latex foley catheter with chemically bonded hydrophilic coating, internally and externally, to latex core.
- Easy and comfortable insertion, hydrophilic surface cushions catheter against urethra.
- Hydrogel resists bacterial colonisation reducing risks of UTI.
- Resists encrustation and colonisation.