


CATHETER COLOPLAST 14 FG MENS INTERMITTENT

\$2.00



SKU: 10250040

PRODUCT DESCRIPTION

When your bladder isn't working as well it could be, your doctor may recommend intermittent self-catheterisation to help you empty your bladder of urine.

For this procedure, you'll need a catheter. The Coloplast Self-Cath is a single-use catheter designed for self-catheterisation. It's made from flexible medical grade PVC with smooth fire polished eyelets for comfort and ease of use. Smooth insertion is also aided by a siliconised surface.

Coloplast Self-Cath features:

- 40cm
- Straight tip
- 14G
- 100% latex free