

## CENTRIUS FIT BARIATRIC MATTRESS SYSTEM



SKU: N/A



## PRODUCT DESCRIPTION

### Purpose & Benefits of Centrius FIT Bariatric Mattress System

The **Centrius FIT Bariatric Mattress System** is a high-capacity dynamic support surface designed specifically for bariatric patients at high to very high risk of pressure injuries. Combining advanced pressure redistribution technology with reinforced construction, the system delivers effective therapy while ensuring comfort, safety, and stability for larger users.

With alternating pressure and continuous low pressure modes, the Centrius FIT provides tailored clinical support, while its intelligent auto-weight sensing technology ensures optimal performance without manual adjustment. Designed for use in hospitals, bariatric care facilities, and complex home care settings, this mattress system supports improved outcomes and caregiver efficiency.

### Key Features

- Bariatric dynamic mattress system with high weight capacity
- Alternating pressure and continuous low pressure therapy modes
- Intelligent auto-weight sensing for responsive pressure adjustment
- Reinforced cell-on-cell design for patient safety and system reliability
- Vapour-permeable, waterproof cover supports infection control

- Quiet operation suitable for long-term and overnight use
- Ideal for high to very high-risk bariatric patients in all care settings

## Specifications

Product Name	Centrius FIT
Product Category	Mattress Replacement
Risk	Very High
Waterlow Scale	20+
Pressure Sore Grading	4
Safe Working Load	250kg
Dynamic Mode	✓
Static Mode - Constant Low Pressure	✓
Intelligent Pressure Sensor	✓
Low Pressure/Power Failure Alarms + Battery Backup	✓
Number of Air Cells	20 x 8" square
CPR Fast Deflation	Dial Turn
Cover - 2 Way Stretch	✓
Dartex Cover - 4 Way Stretch	Optional
Power Cable Management System	✓
Washable Cover	✓
Digital Control	✓
True Low Air Loss	-
Integrated Foam or Air Base	Optional
Cycle Times (mins)	10
Upright Mode	✓
Firm Mode	✓
Memory Recall	✓
Auto Panel Lock	✓
Battery Backup for Alternating Mode	-
Mattress Size - Single	✓
Mattress Size - Euro King Single	-
Mattress Size - King Single	✓
Mattress Size - Extra Wide	-
Mattress Size - Double	-
Mattress Size - Queen	-

Mattress Size - Other

-

Mattress Firm Side Bolsters

-