




MULTIPURPOSE STEP STOOL WITH HAND RAIL



SKU: JAN-JMC-5710

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Our Step Stool is a small, multipurpose stepping stool for all the hard-to-reach locations in your home. It is a perfect tool for children, elderly individuals, the disabled, and anyone else who needs the extra boost. The single-step stool measures 9" in height, and it can be used indoors and outdoors on any level surface.

This stool's solid frame is extremely sturdy, so there is no risk of falling off. Each leg is reinforced with anti-skid rubber to prevent sliding or tipping while in use. The non-slip tread on the top platform is designed to provide additional security. No matter what purpose it is used for, the Clinton Step Stool T-40 is guaranteed to last for years.

No need to piece together dozens of small, individual parts. This step stool is made from a single welded frame, so no screws or rivets are necessary for assembly. It also features a sleek chrome finish that prevents scratches and rusting. This stool is easily cleaned with a damp rag and mild soap, making it an ideal fit for outdoor use.

STEPPING STOOL SPECIFICATIONS

- Height of Step: 23cm
- Overall Width: 44cm
- Overall Depth: 35.5cm
- Step Width: 36cm
- Step Depth: 29cm
- Max User Weight: 160kg

Buy this Step Up Stool online at Medimart with AfterPay for your shopping convenience.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au