

**COBALT HEALTH 2LIFTU LIFT UP DINING CHAIR - 48CM****\$1,314.50****SKU: C48-LIFTU-H****PRODUCT DESCRIPTION****Purpose and benefits**

The 2LiftU Lift-Up Dining Chair is designed to assist individuals with reduced mobility in safely transitioning between sitting and standing positions. Its integrated lift mechanism enhances independence and reduces reliance on carers, making it ideal for home or aged care use. The 48 cm width suits standard-sized users who require discreet assistance in dining settings.

**Features**

- Integrated lift mechanism with six adjustable support levels
- Optional Plus mechanism for smoother lift control
- High-density foam seat for comfort and pressure relief
- Solid timber frame available in multiple wood finishes
- Glides available in nylon, felt, or rubber
- Designed and manufactured in the EU

**Specifications**

- Seat width: 48 cm
- Seat depth: 44 cm
- Seat height: 91 cm
- Safe working load (SWL): 125 kg
- Warranty: 12 months