


**COBALT HEALTH 2LIFTU LIFT UP DINING CHAIR - 58CM****\$1,430.00****SKU: C58-LIFTU-H****PRODUCT DESCRIPTION****Purpose and benefits**

The 58 cm version of the 2LiftU Dining Chair offers a wider seat base while retaining the same reliable lift-assist mechanism. It provides enhanced space and comfort for users requiring broader seating support, without compromising mobility or style.

**Features**

- Six-level lift system for assisted standing
- Optional Plus mechanism for enhanced control
- Supportive high-density foam cushioning
- Sturdy hardwood construction with quality finish options
- Floor-friendly glides: nylon, felt, or rubber
- Built in the EU using premium materials

**Specifications**

- Seat width: 58 cm

- Seat depth: 44 cm
- Seat height: 91 cm
- Safe working load (SWL): 250 kg
- Warranty: 12 months