


**COBALT HEALTH 2LIFTU LIFT UP DINING CHAIR - 68CM****\$1,430.00****SKU: C68-LIFTU-H****PRODUCT DESCRIPTION****Purpose and benefits**

This 68 cm wide version of the 2LiftU chair is ideal for users needing extra space or enhanced lateral support. It combines durable construction with a smart lift mechanism to improve user independence and reduce strain during transfers.

Features

- Lift mechanism with six adjustable levels for standing assistance
- Optional Plus control feature for progressive support
- Contoured seat with high-resilience foam
- Solid wood frame ensures long-term durability
- Glide options suited for all floor types
- Manufactured in Europe to high standards

Specifications

- Seat width: 68 cm

- Seat depth: 44 cm
- Seat height: 91 cm
- Safe working load (SWL): 250 kg
- Warranty: 12 months