



## **COBALT HEALTH BARIATRIC AIR PRESSURE CARE CHAIR**

**\$2,955.00**



**SKU: C01-BARIATRIC-BLU**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Cobalt Health Air Chair Bariatric is purpose-built to deliver exceptional comfort, support, and mobility assistance for larger individuals requiring extended seated care. It features a reinforced frame, a wider seat, and increased structural stability—while still offering the pressure relief and postural support benefits of the standard Air Chair range.

With tilt-in-space and recline functionality, this chair promotes improved circulation and comfort for long-term use. Carers benefit from its smooth-rolling castors and ergonomic adjustments, while users enjoy superior seating support with padded cushions, headrest, and footrest. Ideal for bariatric care in hospitals, aged care facilities, or home environments.

### **Features**

- Reinforced bariatric design for larger users
- Tilt-in-space and reclining seat functionality
- Pressure-relieving foam in seat and backrest
- Wide, high-back design with lateral support
- Swing-away padded footrest with height adjustment

- Padded armrests and removable headrest
- Smooth manoeuvrability with locking castors
- Optional meal tray available
- Suitable for hospital, aged care, and home care use

## **Specifications**

- Safe Working Load: 200 kg
- Seat Height: 700 mm
- Seat Width: 650 mm
- Seat Depth: 560 mm
- Overall Height: 1300 mm
- Overall Width: 850 mm
- Overall Depth: 1250 mm
- Seat Back Adjustment Angle: 97° – 136°
- Footrest Adjustment Angle: 0° – 90°
- Seat Tilt Adjustment Angle: 0° – 17°
- Included Accessories: Headrest, Footrest Pads
- Optional Accessory: Meal Tray